

Note: This chart is both suitable for Filter Pipette Tips and non-filtered tips

Brand Name	Series	Channel	Volume(μL)	10μL	10μL Extended	20μL	100μL	200μL	300μL	1000μL	1250μL	
Axygen	Axyjet	Single Channel	0.1-2	√	√							
			0.5-10	√	√							
			2-20			√	√	√	√			
			5-50			√	√	√	√			
			10-100			√	√	√	√			
			20-200			√	√	√	√			
		100-1000								√	√	
		Multi-Channel	1-10	√	√							
			5-50			√	√	√	√			
			20-200			√	√	√	√			
50-300				√	√	√	√					
Biohit	Proline	Single Channel	0.1-2.5	√	√							
			0.2-10	√	√							
			0.5-10	√	√							
			5-50			√	√	√	√			
			5-100			√	√	√	√			
			10-250			√	√	√	√			
			10-500								√	√
			20-200			√	√	√	√			
			50-1000								√	√
			100-1000								√	√
	Multi-Channel	0.2-10	√	√								
		5-100			√	√	√	√				
		25-250			√	√	√	√				
		50-250			√	√	√	√				
		5-100			√	√	√	√				
	e-Line	Single Channel	0.1-10	√	√							
			0.2-10	√	√							
			5-120			√	√	√	√			
			10-300			√	√	√	√			
			50-1000								√	√
Multi-Channel		0.2-10	√	√								
		5-120			√	√	√	√				
		10-300			√	√	√	√				
m-Line	Single Channel	0.1-3	√	√								
		0.5-10	√	√								
		10-100			√	√	√	√				
		20-100			√	√	√	√				
		100-1000								√	√	
	Multi-Channel	0.5-10	√	√								
		5-100			√	√	√	√				
		30-300			√	√	√	√				

Thermo Scientific	Finnpipette	Single Channel	0.2-2	√	√						
			0.3-3	√	√						
			0.2-10	√	√						
			0.5-5	√	√						
			0.5-10	√	√						
			1-10	√	√						
			5-40			√	√	√	√		
			5-50			√	√	√	√		
			10-100			√	√	√	√		
			40-200			√	√	√	√		
			30-300			√	√	√	√		
			100-1000							√	√
			200-1000							√	√
	Multi-Channel	1-10	√	√							
		5-50			√	√	√	√			
		50-300			√	√	√	√			
Finnpipette Digital	Single Channel	40-200			√	√	√	√			
Finnpipette electronic	Multi-Channel	1-10	√	√							
		30-300			√	√	√	√			
		100-1000							√	√	
Gilson Pipetman	Pipetman Classic	Single Channel	P2,0,1-2	√	√						
			P2,0,2-2	√	√						
			P10,0,2-10	√	√						
			P20,2-20			√	√	√	√		
			P100,10-100			√	√	√	√		
			P200,20-200			√	√	√	√		
			P1000,100-1000							√	√
	Pipetman Neo	Single Channel	P2N,0,2-2	√	√						
			P10N,1-10	√	√						
			P20N,2-20			√	√	√	√		
			P100N,10-100			√	√	√	√		
			P100N,20-100			√	√	√	√		
			P200N,20-200			√	√	√	√		
		Multi-Channel	P8×20N,2-20			√	√	√	√		
	Pipetman Ultra	Single Channel	U200,20-200			√	√	√	√		
		Multi-Channel	P200-M8,20-200			√	√	√	√		